Warum testen?

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
Niedriges Risiko

PROLARIS MOLECULAR SCORE
A measure of cell proliferation, independent of clinical features.

1 2.6

Prolaris Molecular Score is calculated by directly measuring the biology of the tumor through the expression of cell cycle progression genes.
Based on analysis from two cohorts of conservatively managed patients,^{1,2} Prolaris Molecular Scores between 1.9 and 8.7 are clinically validated and will be reported as the calculated Prolaris Molecular Score.
Prolaris Molecular Score is combined with the patient's clinical features in a clinically validated weighted algorithm resulting in patient's 10-year prostate cancer disease specific mortality (DSM) and patient's 10-year metastasis (METS) risk.

Disease Specific Mortality (DSM) Risk is within the threshold for active surveillance*

FEATURES USED FOR RISK ASSESSMENT

Prolaris Molecular Score: 2.6
Patient Age at Biopsy: 66
PSA prior to this Biopsy: 7.1 ng/ml **2**
Clinical T Stage: T1c
% Positive Cores: 33%
Gleason Score: 3+4=7 (Group 2 ISUP³)
CAPRA Score⁴: 3

RISK CLASSIFICATION BASED ON GUIDELINES

NCCN Risk Category⁵: Favorable intermediate

Hohes Risiko

PROLARIS MOLECULAR SCORE
A measure of cell proliferation, independent of clinical features.

1 5.2

Prolaris Molecular Score is calculated by directly measuring the biology of the tumor through the expression of cell cycle progression genes.
Based on analysis from two cohorts of conservatively managed patients,^{1,2} Prolaris Molecular Scores between 1.9 and 8.7 are clinically validated and will be reported as the calculated Prolaris Molecular Score.
Prolaris Molecular Score is combined with the patient's clinical features in a clinically validated weighted algorithm resulting in patient's 10-year prostate cancer disease specific mortality (DSM) and patient's 10-year metastasis (METS) risk.

Disease Specific Mortality (DSM) Risk exceeds the threshold for active surveillance*

FEATURES USED FOR RISK ASSESSMENT



Prolaris Molecular Score: 5.2
Patient Age at Biopsy: 67
PSA prior to this Biopsy: 7.2 ng/ml **2**
Clinical T Stage: T1c
% Positive Cores: 33%
Gleason Score: 3+4=7 (Group 2 ISUP³)
CAPRA Score⁴: 3

RISK CLASSIFICATION BASED ON GUIDELINES

NCCN Risk Category⁵: Favorable intermediate

Prolaris Molecular Score and clinical features are combined in a clinically validated weighted algorithm

Prolaris Molecular Score and clinical features are combined in a clinically validated weighted algorithm

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1. Eppstein JI, et al. Am J Surg Pathol 2014 2. NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines®) for Prostate Cancer v2.2020. © National Comprehensive Cancer Network, Inc. 2020. All rights reserved. Accessed November 10, 2020 2. 3. Cuzick J, et al. Lancet Oncol 2011 4. Cuzick J, et al. BJC 2012 5. Cooperberg MR, et al. J Clin Oncol 2013 6. Freedland SJ, et al. Int J Radiat Oncol Biol Phys 2013 7. Bishoff JT, et al. J Urol 2014 8. Cuzick J, et al. Br J Cancer 2015 9. Koch MO, et al. Cancer Biomarkers 2016 10. Tosoian J, et al. BJU International 2017 11. Lin DW, et al. Urologic Oncology 2018 12. Canter DJ, et al. European Urology 2018 13. Leon P, et al. World Journal of Urology 2018 14. Canter DJ, et al. Prostate Cancer and Prostatic Diseases 2019 15. Tward et al. Clinical Genitourinary Cancer 2021 16. et al. The Prostate 2021